


RH01687

Product Code: F311500
 CAS No.: 302901-13-9
 Mol. Weight: C12H9ClN6O2S
 Mol. Formula: 336.76
 Storage Temp: Powder: -20°C, 36 months; 4°C, 24 months
 In solvent: -80°C, 6 months; -20°C, 1 month
 Target: Others
 Pathway: Others



Biological Activity

制备储备液:

质量 溶剂体积 浓度	1 mg	5 mg	10 mg
1 mM	2.9695 mL	14.8474 mL	29.6947 mL
5 mM	0.5939 mL	2.9695 mL	5.9389 mL
10 mM	0.2969 mL	1.4847 mL	2.9695 mL

Bioactive: RH01687 is a compound that can protect pancreatic β cells against endoplasmic reticulum stress-induced cell death. RH01687 has the potential for the research of diabetes.



访问 FELIX 中国官网获取最新优惠活动、产品动态
FELIX 所有产品只能用于科学研究
任何个人和单位不能将产品用于人体、临床等非科学研究领域和用途