

## Glucocorticoids receptor agonist 3

Product Code: F312890

CAS No.: 305821-96-9


Mol. Weight: C31H32N2O2

Mol. Formula: 464.6

Storage Temp: Please store the product under the recommended conditions in the Certificate of Analysis.

Target: Glucocorticoid Receptor

Pathway: Immunology/Inflammation; Vitamin D Related/Nuclear Receptor



## Biological Activity

制备储备液:

质量 溶剂体积 浓度	1 mg	5 mg	10 mg
1 mM	2.1524 mL	10.7619 mL	21.5239 mL
5 mM	0.4305 mL	2.1524 mL	4.3048 mL
10 mM	0.2152 mL	1.0762 mL	2.1524 mL

**Bioactive:** Glucocorticoids receptor agonist 3 is a potent agonist of glucocorticoids receptor. Glucocorticoids receptor agonist 3 is useful to research diseases, such as obesity, diabetes, and inflammation (extracted from patent W02000066522A1, compound 345).



访问 FELIX 中国官网获取最新优惠活动、产品动态  
FELIX 所有产品只能用于科学研究  
任何个人和单位不能将产品用于人体、临床等非科学研究领域和用途